

Med-O-Wheel GSM as a part of a IoT -service

Med-O-Wheel GSM works as the Smart, but also allows to follow the medication online. You can see the taken or not taken medicine remotely.

The Addoz Med-O-Wheel Smart, when connected to a IoT - system, is not geographically restricted and safe medication can be carried out anywhere.

IoT- function requires Med-O-Wheel Smart equipped with GSM-module.

Addoz Med-O-Wheel Smart

- can be timed to dispense 1-4 doses per day
- lockable
- using dosage cassettes with 14 and 28 dosage compartments
- weight 296 g
- dimensions \varnothing 120 x 132 x 60 mm
- battery type 2 x AA alkaline batteries
- battery life about one year



Addoz Med-O-Wheel

- easy-to-use medication dispenser
- available in 14 and 28 dosage compartments
- shock-resistant and machine washable
- serves the doses one compartment at a time
- cover available in a range of colours
- weight 110 g
- dimensions \varnothing 125 x 35 mm



AZ10127-0



Addoz Med-O-Wheel Smart - the intelligent medication dispenser



“Have I already taken my medicine?”

is a question that daily occupies the minds of those who require medicine. The worry of remembering to take the medicine is shared by close relatives and caregivers. According to expert estimates, only one-third of medication is taken as instructed, one-third randomly and one-third not at all.

The neglect of long-term medication causes a significant risk in institutional care. Through the safe and correctly timed medicine administration by the Addoz Med-O-Wheel Smart, you can enhance your well-being and vitality and improve your chance of living at home for a longer period of time.



Addoz Med-O-Wheel Smart

The Addoz Med-O-Wheel Smart is a new, easy-to-use dispenser for medicines in pill and capsule form. The Med-O-Wheel Smart will help you to administer your medicine at the right time, and prevent unintentional misuse. The device ensures that the medicine is taken at the right time, and it reminds by audible and light signals when the user ignores to take the medicine. When the user presses the transparent cover of the Med-O-Wheel Smart, the device rotates the correct dose into the accessible position, and then waits for the next specified dosage time.

If, despite everything, a dose is missed, the Med-O-Wheel Smart will move the skipped dose aside for monitoring. The device only gives the new dose at the right time, which in practice prevents the deliberate or accidental misuse of the medicine.

You can easily put the medicine into the dosage cassette yourself, but this can also be done at the chemist's. The chemist will then send the cassette to the user or care institution. An emptied dosage cassette is correspondingly returned to the chemist for refilling. The cassette is very suitable for medication logistics.

Features and benefits

Doses of medicine only available at the right time
→ *improves user safety*

No power cords or wires
→ *can be taken along and used anywhere*

Dosage cassette easy to fill
→ *can be done by the user, close relative or professional carer*

Alarm sounds if a dose is missed
→ *reduces the problems caused by inexact administering of medicine*

Automatic low battery warning
→ *ensures continuous optimal operation*

In order to take the medicine, the user must press the button on the dispenser
→ *monitors the activity of the user*

Suitable for all

The Addoz Med-O-Wheel Smart is suitable for people of all ages, from families with children to senior citizens. It relieves the user and his/her close relatives from the worry of remembering to administer the medicine. It ensures that the medicine is taken as instructed, both in long and short term courses of medication. The Med-O-Wheel Smart is also suitable for use by a care organisations GSM-enabled remote care service.